

# TRANSPARENT MASTIC

MARBLE • GRANITE  
ONYX • STONE



## TECHNICAL DATA SHEET

### DESCRIPTION

Mastic for marble, granite, onyx and stone. Once hardened, the product is dry and well polishable. It has good storage stability.

Can be easily coloured by addition of paste pigment for polyesters or metal oxides.

### SAFETY RULES

Information provided in this technical data sheet is based upon our experience. We assure the perfect quality of products. However, being the utilization out of our control, we don't assume responsibility on the result obtained. It is obligation of the customer to verify that the use of the product is done according to laws, legislation, patents of the customer's country.

### INSTRUCTIONS OF USE

- Before use the glue check that the surface has to be dry, clean and dust free.
- Use a spatula for take the necessary quantity of glue from the can and add 1-3% of the hardener in paste
- Pay attention to not add the hardener in paste inside the can glue.
- Mix well and use the product. Do not return the unused glue in the can
- For colour the glue use Uni Colorant Paste, add the colour to the glue, max 3%, before and then add the hardener in paste.

### SHELF LIFE

The products will last at least 12 months in the normal condition between 18-25°C, kept away from source of heat, humidity and sun light.

# TRANSPARENT MASTIC

MARBLE • GRANITE  
ONYX • STONE



## TECHNICAL DATA

	Solid	Semisolid	Liquid
Specif weight	1,78 g/cm <sup>3</sup> ±0,05	1,75 g/cm <sup>3</sup> ±0,05	1,12 g/cm <sup>3</sup>
V.O.C Directive 2010/75/EC	12,68% - 228,24 gr/litre	12,68% - 228,24 gr/litre	12,68% - 228,24 gr/litre
Viscosity cps 25°	Tixotropic paste	200.000 - 300.000	1300 - 1600
Consistency/Appereance	Vertical colored paste	Cream colored paste	Liquid transparent paste
Application	Reparing/Bonding	Filling big holes	Filling holes and cracks
Hardener	1-3%	1-3%	1-3%
Dry to touch at 25°	15 min	15 min	15 min
Gel time at 25°	3-4 min	3-4 min	5-8 min
Complete curing at 25°	20 min	20 min	20 min
Polishable at 25°	30 min	30 min	30 min